

Active®



SPIRIT

Laser Hair Removal

Beauty Starts Here




[www.active-lasers.com](http://www.active-lasers.com)



**High Power and Speed – The innovative hair removal solution for all skin tones, sensitive and aged skin.**

Active Lasers **SPT** – (Speed Pulse Technology) is the winning acquisition for aesthetics centers around the world.

Spirit's **SPT** technology uses high quality laser diodes to deliver high intensity short pulses for noticeable results from the first treatment.



# Spirit

The perfect solution for laser hair removal



**Cold Touch Hand Piece**  
Continuous contact cooling for virtually pain-free treatments.



**Fast Treatments**  
With up to 10 pulses per second, treatment time is significantly reduced.



**Superior Safety Standards**  
Spirit's hand piece is activated only while in direct skin contact.



**Three Wavelengths**  
755nm, 810nm and 1064nm enabling treatment for all skin tones and hair colors.



**Fast ROI**  
Exceptional technology and diode precision offset the cost via rapid revenue



**Cooling System**  
A powerful cooling system allows operation at any room temperature.



**Outstanding results From the first treatment**

## Technical Specification

<b>Cooling System</b>	Maintains optimum temperature regardless of room temperature
<b>Hand Piece</b>	Provides a cooling sensation for a virtually painless treatment
<b>Pulses per second</b>	1-10 Pulses Per Second
<b>Combined Wavelengths</b>	755nm, 810nm, 1064nm
<b>Output Power</b>	3000 Watt
<b>Treatment area</b>	(20x14mm) 2.8cm <sup>2</sup>
<b>Built in safety</b>	Hand piece activated by direct skin touch
<b>Weight</b>	42 Kg
<b>Size</b>	44cm x 45cm x 128cm
<b>Interface</b>	Computer touch screen

## Hand Piece

The ergonomic design enables easy hair removal including in difficult to reach, intimate and sensitive areas.

Ability to customize treatments for superb results. The SPT laser targets the dermal papilla, inhibiting hair growth from the very first treatment.

Spirit's hand piece is equipped with a built-in safety mechanism activating only while the hand piece is in direct skin contact.



Before Treatment



After Treatment





# WHY ACTIVE?

Our innovative products will help you grow your customer's base by providing solutions in line with the needs of the modern aesthetic clients.



## BECOME A PART OF SUCCESS

Expand your portfolio with Advanced Medical Technology

- +972-3-923-6610
- marketing@active-lasers.com
- www.active-lasers.com
- Shoham 11 st, Petah Tikva, Israel. 4951766



Global Positioning Across  
**52**  
Countries and regions  
Worldwide



More Than  
**5,926**  
Devices Installed  
Worldwide



Over  
**40**  
Years of Experience



Attractive price with low  
maintenance and service  
cost.

## Beauty Starts Here